BOOST 2 CORE

effortless off-weighting.





KEY FEATURES

Percentage of Body Weight	100% - 20%
Incremental Adjustment	1% Increments
Height Range	4'6" to 7' (137 to 208 cm)
Clinical Gait Views	Front, Back, Side
23" Touchscreen with Knob Control	•
Pre-programmed Workouts	•
Shorts (XS-XXL)	One Pair Included
10 Pack Shorts	Optional
Speed	12 MPH (19 KM/h)
Incline	15%
Length	82" (208 cm)
Width	42" (107 cm)
Color	White
115 Vac 20 A	•
208/230 Vac 20 A	Optional
Commercial Warranty	1 Years Parts and Labor

Contact us at 800-WOODWAY to speak to a team member for additional information.

All specifications contained herein is based on info at time of publication and subject to change.

Rev 230711

SIMPLIFIED USER CONTROL

No longer do you need to press and hold buttons - the Boost 2 offers quick and accurate adjustment of speed and body weight via an easy and intutive knob control.

PRECISION HEIGHT CONTROL

Precision height control lifts precisely from 4'6" to 7'.

With no manual lifting involved, the Boost 2 offers a quick, easy, and safe, entry/exist.

SLAT BELT

The original patented Woodway slat-belt remains best-in-class. The Woodway belt is clinically proven to reduce shock and increase comfort.





