

## BARIATRIC SOLUTION

### THERAPY TREADMILL TO ACCOMMODATE A VARIETY OF PATIENTS

Much like the PPS Med, the Bari-Mill provides comfort and assistance to patients who require walking support. With infinite possibilities in width and height the Bari-Mill is great for larger bariatric patients as well as children due to the lowering capabilities.



*\*Various dimensions available*

For The Long Run®

**WOODWAY®**

WOODWAY WORLD HEADQUARTERS  
W229 N591 Foster Court Waukesha, WI 53186  
PHONE (800) 966 - 3929 | FAX (262) 522 - 6235  
**WOODWAY.COM**

WOODWAY INTERNATIONAL GMBH  
Steinackerstrasse 20, D-79576 Weil am Rhein, Germany  
PHONE +49 - 7621 - 940 999 - 0 | FAX +49 - 7621 - 940 999 - 40  
**WOODWAY.DE**

## BARI-MILL

- Standard Feature  
○ Optional Feature



USER INTERFACE	BARI-MILL	
	Coded/Non-Coded Polar Heart Rate Pick-up	●
	ANT+ Heart Rate Pick-up	
	USB Power Supply	
	Touch Grip Heart Rate Pick-up	
	Woodway User System (WUS)	●
	RS 232 Serial Interface	●
	CSAFE Communication Compatible	
	Wired Remote	●
	Running Surface Dimensions	22" W x 68" L (55 x 173 cm)
SPECIFICATIONS	User Weight Capacity	Run: 400 lb. / Walk: 800 lb. (4 mph max)
	Belt Type	60 Individual Slats
	Drive System	114 Precision Ball Bearings with 12 Roller Guides (4 mm lateral tolerance)
	Running Surface	Vulcanized Rubber (38-43 shore hardness)
	Drive Motor	2 hp Continuous (5 hp peak) Brushless Servo
	Unit Weight	490 lb. (Shipping Weight 575 lb.)
	Width <i>*Alternative Treadmill Frame Widths Options Available</i>	41" (104 cm)
	Length	73" (185 cm)
	Height	62" (158 cm)
	0-25% Incline	●
	(-5%) - (+20%) Incline	○
	0-12.5 MPH (0-20 km/h)	●
	0-15 MPH (0-24 km/h)	○
	Reverse	●
	115 Vac 20 Amp	●
	208/230 Vac 20 Amp	○
	Hospital Grade Circuitry	●
	Therapist Side Seating	○
	Warranty	5 Year Drive, Motor, & Belt 4 Year All Components 1 Year Labor

*Treadmill specifications subject to change without written notice.  
Warranties outside the U.S. may vary - Please contact your dealer for details.*

*Full Compliance with the Legal Requirements for  
Medical Products in the European Economic Area (CE marking)*

*WOODWAY PPS Series treadmills comply with the following:*

- Medical Devices Directive (directive 93/42/EEC)\*
- Medical Standard EN 60601-1, 3rd edition
- Medical Standard UL 60601-1 (USA)
- Medical Standard UL CSA C22.2 No 601.1 (Canada)



## CONVENIENT USER CONSOLE

- LCD screen with readouts monitoring speed, incline, distance, time, and heart rate
- Hand-held remote (speed, elevation & stop)



## WOODWAY SLAT BELT RUNNING SURFACE

- Individual slats with vulcanized rubber surface
- Zero stretch & Zero slip toothed, wire-reinforced lateral belts
- Near frictionless precision ball bearing rail guide

## SAFETY &amp; RELIABILITY

- Emergency stop magnet and safety lanyard
- Hospital grade circuitry
- Polar monitor circuitry – integrated heart rate monitoring
- Advanced Slat Belt shock absorption
- Life expectancy of running surface 150,000 to 200,000 miles
- Extended warranty and factory renewal program

**FOR ADDITIONAL INFORMATION ABOUT THIS PRODUCT  
VISIT [WWW.WOODWAY.COM](http://WWW.WOODWAY.COM) OR CALL 800.WOODWAY TO TALK TO A REPRESENTATIVE.**