

HEALING ONE STEP AT A TIME

REHABILITATION TREADMILL

From rehabilitation to orthopedics and physical therapy, the PPS Med (Professional Power System) treadmill provides comfort and assistance to patients who require walking support. With an absolute zero start speed in increments of 0.1 mph a patient can feel assured that they are in control and can set their own pace that best suits their comfort levels.



**Various widths available*

For The Long Run®

WOODWAY®

WOODWAY WORLD HEADQUARTERS
W229 N591 Foster Court Waukesha, WI 53186
PHONE (800) 966 - 3929 | FAX (262) 522 - 6235
WOODWAY.COM

WOODWAY INTERNATIONAL GMBH
Steinackerstrasse 20, D-79576 Weil am Rhein, Germany
PHONE +49 - 7621 - 940 999 - 0 | FAX +49 - 7621 - 940 999 - 40
WOODWAY.DE

PPS MED



- Standard Feature
- Optional Feature

PPS MED		
USER INTERFACE	Coded/Non-Coded Polar Heart Rate Pick-up	●
	ANT+ Heart Rate Pick-up	
	USB Power Supply	
	Touch Grip Heart Rate Pick-up	
	Woodway User System (WUS)	●
	RS 232 Serial Interface	●
	CSAFE Communication Compatible	
	Wired Remote	●
	Running Surface Dimensions	22" W x 68" L (55 x 173 cm)
	User Weight Capacity	Run: 400 lb. / Walk: 800 lb. (4 mph max)
SPECIFICATIONS	Belt Type	60 Individual Slats
	Drive System	114 Precision Ball Bearings with 12 Roller Guides (4 mm lateral tolerance)
	Running Surface	Vulcanized Rubber (38-43 shore hardness)
	Drive Motor	2 hp Continuous (5 hp peak) Brushless Servo
	Unit Weight	465 lb. (Shipping Weight 533 lb.)
	Width	41" (104 cm)
	*Alternative Treadmill Frame Widths Options Available.	
	Length	73" (185 cm)
	Height	59" (150 cm)
	0-25% Incline	●
	(-5%) - (+20%) Incline	○
	0-12.5 MPH (0-20 km/h)	●
	0-15 MPH (0-24 km/h)	○
	Reverse	●
	115 Vac 20 Amp	●
	208/230 Vac 20 Amp	○
	Hospital Grade Circuitry	●
	Therapist Side Seating	○
	Parallel Bars	○
	Front Handrail	Call
Warranty	5 Year Drive, Motor, & Belt 4 Year All Components 1 Year Labor	



CONVENIENT USER CONSOLE

- LCD screen with readouts monitoring speed, incline, distance, time, and heart rate
- Hand-held remote (speed, elevation & stop)



WOODWAY SLAT BELT RUNNING SURFACE

- Individual slats with vulcanized rubber surface
- Zero stretch & Zero slip toothed, wire-reinforced lateral belts
- Near frictionless precision ball bearing rail guide

SAFETY & RELIABILITY

- Emergency stop magnet and safety lanyard
- Hospital grade circuitry
- Polar monitor circuitry – integrated heart rate monitoring
- Advanced Slat Belt shock absorption
- Life expectancy of running surface 150,000 to 200,000 miles
- Extended warranty and factory renewal program

Treadmill specifications subject to change without written notice. Warranties outside the U.S. may vary - Please contact your dealer for details.

Full Compliance with the Legal Requirements for Medical Products in the European Economic Area (CE marking)

- WOODWAY PPS Series treadmills comply with the following:
- Medical Devices Directive (directive 93/42/EEC)*
 - Medical Standard EN 60601-1, 3rd edition
 - Medical Standard UL 60601-1 (USA)
 - Medical Standard UL CSA C22.2 No 601.1 (Canada)

**FOR ADDITIONAL INFORMATION ABOUT THIS PRODUCT
VISIT WWW.WOODWAY.COM OR CALL 800.WOODWAY TO TALK TO A REPRESENTATIVE.**